



## THE WORLD LEADER IN DISPOSABLE DIGESTION TUBES

RackLock FEATURE ALLOWS ONE-HANDED CLOSURE OF SCREW CAP

DISPOSABLE - ELIMINATES THE NEED TO CLEAN BEAKERS

CLASS A 50 ML & 100 ML GRADUATIONS FOR EASY & ACCURATE DILUTIONS

ULTRA LOW LEACHABLE METAL CONTENT TO ELIMINATE CONTAMINATION

LEAK-FREE SCREW CAP

USED IN DIGIPREP OR OTHER GRAPHITE BLOCK DIGESTION SYSTEMS



This single use, all in one, universal vessel saves time, money and eliminates the need to stock 5 different vessels!

### The revolutionary *DigiTUBE* replaces:

- Volumetric flasks
- Graduated cylinders
- Autosampler tubes
- Storage containers
- Digestion vessels

### Did you know about the *DigiTUBE* Calibration Certificate?

Every case of *DigiTUBE*s comes complete with a lot specific Calibration Certificate that guarantees Class A tolerance at the 50 ml and 100 ml\* line of every *DigiTUBE*. Quality, consistency and cost effectiveness... that's the **SCP SCIENCE** way.

\* for 100 ml *DigiTUBE*.

### Also Available for the Cleanroom

*DigiTUBE*s are available in shrink wrap packaging to eliminate contamination from storage and handling.

## ORDERING INFORMATION

	Catalogue #	Description	Quantity / Packaging
Bulk <i>Digi</i> TUBEs	010-515-261	<i>Digi</i> TUBEs 15 ml RL w/caps	1620
	010-515-262	<i>Digi</i> TUBEs 15 ml RL w/o caps	1620
	010-515-263	<i>Digi</i> TUBEs 15 ml NRLock w/caps	1620
	010-515-264	<i>Digi</i> TUBEs 15 ml NRLock w/o caps	1620
	010-515-361	<i>Digi</i> TUBEs 15 ml RL w/caps	540
	010-515-363	<i>Digi</i> TUBEs 15 ml NRLock w/caps	540
	010-500-261	50 ml <i>Digi</i> TUBE with RackLock, blue caps	750
	010-500-262	50 ml <i>Digi</i> TUBE with RackLock, no caps	750
	010-500-263	50 ml <i>Digi</i> TUBE without RackLock, blue caps	750
	010-500-264	50 ml <i>Digi</i> TUBE without RackLock, no caps	750
	010-501-261	100 ml <i>Digi</i> TUBE with RackLock, blue caps	300
	010-501-262	100 ml <i>Digi</i> TUBE with RackLock, no caps	300
	010-501-263	100 ml <i>Digi</i> TUBE without RackLock, blue caps	300
	010-501-264	100 ml <i>Digi</i> TUBE without RackLock, no caps	300
The 50 ml <b>HandyPAK</b> has 25 <i>Digi</i> -TUBEs shrink wrapped to eliminate contamination during storage or transport (30 packs, blue caps supplied in separate bag). <b>HandyPAK</b> Plus has 24 capped <i>Digi</i> TUBEs in storage racks, shrink wrapped (20 packs)	010-500-362	<b>HandyPAK</b> , 50 ml <i>Digi</i> TUBE with RackLock, blue caps, shrink wrapped	30 x 25
	010-500-364	<b>HandyPAK</b> , 50 ml <i>Digi</i> TUBE without RackLock, blue caps, shrink wrapped	30 x 25
	010-500-461	<b>HandyPAK</b> Plus, 50 ml <i>Digi</i> TUBE with RackLock, blue caps, shrink wrapped	20 x 24
	010-500-463	<b>HandyPAK</b> Plus, 50 ml <i>Digi</i> TUBE without RackLock, blue caps, shrink wrapped	20 x 24
The 100 ml <b>HandyPAK</b> has 9 <i>Digi</i> -TUBEs shrink wrapped to eliminate contamination during storage or transport (34 packs, blue caps supplied in separate bag). <b>HandyPAK</b> Plus has 12 capped <i>Digi</i> TUBEs in storage racks, shrink wrapped (24 packs)	010-501-362	<b>HandyPAK</b> , 100 ml, <i>Digi</i> TUBE with RackLock	9 x 34
	010-501-364	<b>HandyPAK</b> , 100 ml, <i>Digi</i> TUBE without RackLock	9 X 34
	010-501-462	<b>HandyPAK</b> Plus, 100 ml, <i>Digi</i> TUBE with RackLock	24 x 12
	010-501-464	<b>HandyPAK</b> Plus, 100 ml, <i>Digi</i> TUBE without RackLock	24 x 12
Digestion Rack for 15 ml <i>Digi</i> TUBEs	010-510-023	25 position RackLock Digestion Rack for 15 ml <i>Digi</i> TUBE, for use with <i>Digi</i> PREP MINI	1
	010-510-027		1
Digestion Rack for 50 ml <i>Digi</i> TUBEs	010-500-021	24 position RackLock Digestion Rack for 50 ml <i>Digi</i> TUBE, for use with <i>Digi</i> PREP MS/LS	1
	010-500-025	18 position RackLock Digestion Rack for 50 ml <i>Digi</i> TUBE, for use with HotBlock™	1
	010-505-021	24 position RackLock Digestion Rack for 50 ml <i>Digi</i> TUBE, for use with <i>Digi</i> PREP Jr.	1
Digestion Rack for 100 ml <i>Digi</i> TUBEs	010-501-510	12 position RackLock Digestion Rack for 100 ml <i>Digi</i> TUBE, for use with <i>Digi</i> PREP Jr.	1
	010-501-530	15 position RackLock Digestion Rack for 100 ml <i>Digi</i> TUBE, for use with <i>Digi</i> PREP MS	1
	010-501-550	21 position RackLock Digestion Rack for 100 ml <i>Digi</i> TUBE, for use with <i>Digi</i> PREP LS	1
Watch Glasses	010-500-081	Disposable Polypropylene Watch Glasses, 50 ml	1000
	010-501-081	Disposable Polypropylene Watch Glasses, 100 ml	500
Storage Racks	010-500-026	Plastic Storage Racks for 50 ml <i>Digi</i> TUBE	5
	010-501-028	Foam Storage Racks for 100 ml <i>Digi</i> TUBE	5
Screw Caps	010-515-060	Screw Cap, Blue 15 ml <i>Digi</i> TUBE	540
	010-515-172	Screw Cap, Clear, 15 ml <i>Digi</i> TUBE	540
Color Coded 50 ml <i>Digi</i> TUBEs - Choose 3 colors of your choice per <i>Digi</i> TUBE box	010-500-060	Blue Caps for 50 ml <i>Digi</i> TUBEs	250
	010-500-130	Yellow Caps for 50 ml <i>Digi</i> TUBEs	250
	010-500-140	Orange Caps for 50 ml <i>Digi</i> TUBEs	250
	010-500-150	Red Caps for 50 ml <i>Digi</i> TUBEs	250
	010-500-172	Clear Caps for 50 ml <i>Digi</i> TUBEs	250
Color Coded 100 ml <i>Digi</i> TUBEs	010-501-060	Blue Caps for 100 ml <i>Digi</i> TUBEs	100

### Authorised Distributor:

Instrument Solutions Benelux BV  
 De Liesbosch 50 - 3439 LC Nieuwegein (NLD)  
 Tel: +31(0)88 4678 786 / Fax: +31(0)88 4678799  
[www.instrument-solutions.com](http://www.instrument-solutions.com)  
[info@instrument-solutions.com](mailto:info@instrument-solutions.com)